

Extraction Set  
Hip Revision  
without Bone Loss

**GOMINA**  
SWISS QUALITY





## A World First from Gomina

For the first time, Gomina has developed an extraction set that allows a hip revision to be done without bone loss. Thanks to innovative production technologies, it is now possible to manufacture dimensionally rigid chisels with different radii. The carefully grinded, sharp cutting edges of the chisels are the result of Gomina's many decades of expertise and make a precise separation of bone and prosthesis possible.

The newly developed handle with a quick clamping mechanism offers various fixation positions for the chisels and impresses with its ergonomic design. Chisels and handle are perfectly matched, thus ensuring maximum transmission of force during use. This innovative, pioneering extraction set from Gomina makes surgeons' work easier and leads to shorter operating times. Entirely for the patient's benefit.

**GOMINA AG** Raiftstrasse 4, CH-3989 Niederwald

Telefon ++41 (0)27 970 12 70, Fax ++41 (0)27 970 12 75. [www.gomina.ch](http://www.gomina.ch), [info@gomina.ch](mailto:info@gomina.ch)