

**Surgical Saw Blades for Orthopaedics,
Trauma and Small Bones
fitting Aesculap®
PRODUCT CATALOGUE**



REUSABLE
SAW BLADES
(except Keel Blade)



Overview Gomina Saw Blades

Saw blades for Orthopaedics



Saw blades for Trauma



Saw blades for Small Bones



Reciprocating Saw blades



Sternum Saw blades



Reciprocating Keel Blade



Keel Blade Set



Saw Blade Portfolio

Page

Fits to Aesculap®



Aesculap®
Acculan 3TI
GA673, GB660R,
Acculan 4
GA331, GA336



Aesculap®
Acculan 3TI
GB435, GB440,
GA635, GA636,
GA641, GA674
Acculan 4
GA334




Aesculap®
Elan 4 System
GA673, GB660R,
GB891R

10


Teeth designs

Premium teeth design

The Gomina Premium Line of saw blades is characterised by a combination of cross-grinded cutting teeth and straight-grinded raker teeth. This results in optimal cutting efficiency and high precision. The teeth are aligned on a straight line and are asymmetrical, which leads to a high degree of self-centring of the saw blade.



Premium teeth design



Performance teeth design

The Gomina Performance Line of saw blades is characterised by asymmetrical straight-grinded cutting teeth, slightly less aggressive than the Premium Line.




Performance teeth design



Reciprocating teeth designs

The Gomina reciprocating saw blades offers a variety of teeth designs. The cross-grinded cutting teeth (reciprocating aggressive) results in optimal cutting efficiency and high precision.

Reciprocating saw blades with the straight-grinded and symmetrical cutting teeth design are less aggressive. The set teeth allow an excellent sawing performance without jamming. Gomina double-sided reciprocating blades completes the portfolio.




Trauma teeth design

The Gomina trauma saw blades offers the wide product range for the corresponding application. The teeth set allow an excellent sawing performance without jamming.






Small Bone teeth design

The Gomina small bone saw blade is characterised by asymmetrical straight-grinded cutting teeth. These saw blades are mostly used without cutting blocks due the ease for free-hand sawing.



General Saw Blades information



NOTE: All markings on the products in this catalogue are exemplary, positions may differ from the actual marking. The UDI is to be understood as an example. Every article has its individual code.

Packaging

NON-STERILE PACKAGING

Non-sterile saw blades are packaged in single pouches and carry the product label with detailed information (label printed black & white). They don't come in boxes and can be **ordered single** for any desired quantity. Our unsterile saw blades are all reusable.



STERILE PACKAGING

Sterile saw blades are packaged in double sterile pouches which feature an additional protection layer for the cutting edge. Minimal ordering quantity is **one box containing 5 saw blades**. The box also contains one patient label per saw blade, placed between the first and second pouch of each blade. Except for the Reciprocating Keel Blade, all our sterile saw blades are reusable.



Contents

Fits to Aesculap®	9
Premium Fits to: Acculan 3Ti	9
Performance Fits to: Aesculap®	16
Performance Fits to: Acculan 3Ti	17
Reciprocating Fits to: Acculan 3Ti	21
Keel Blade Set	24
Trauma Fits to: Aesculap®	26
Trauma Fits to: Acculan 3Ti	27
Small Bone Fits to: Acculan 3Ti	31
Small Bone Fits to: Elan 4, Acculan 4	34

Fits to Aesculap®

Overview references

Length	Width	Thickness						
Premium Blades		0.90mm	1.00mm	1.18mm	1.19mm	1.23mm	1.27mm	1.37mm
100mm	19mm				219.190.119 219.190.119S			
90mm	25mm	219.259.090 219.259.090S	219.259.100 219.259.100S	219.259.118 219.259.118S	219.259.119 219.259.119S		219.259.127 219.259.127S	219.259.137 219.259.137S
90mm	19/22mm				219.969.119 219.969.119S		219.969.127 219.969.127S	
90mm	19mm	219.199.090 219.199.090S			219.199.119 219.199.119S	219.199.123 n/a	219.199.127 219.199.127S	219.199.137 219.199.137S
90mm	12.5/19mm				219.999.119 219.999.119S		219.999.127 219.999.127S	
90mm	12.5mm	219.129.090 219.129.090S	219.129.100 219.129.100S	219.129.118 219.129.118S	219.129.119 219.129.119S	219.129.123 219.129.123S	219.129.127 219.129.127S	219.129.137 219.129.137S
90mm	6/12.5mm						n/a 219.989.127S	
70mm	12.5mm						219.127.127 n/a	

Note: sterile version with S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.19mm	219.190.119	219.190.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	0.90mm	219.259.090	219.259.090S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	219.259.100	219.259.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.18mm	219.259.118	219.259.118S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.19mm	219.259.119	219.259.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	219.259.127	219.259.127S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.37mm	219.259.137	219.259.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19/22mm	1.19mm	219.969.119	219.969.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19/22mm	1.27mm	219.969.127	219.969.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	0.90mm	219.199.090	219.199.090S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	219.199.119	219.199.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.23mm	219.199.123	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	219.199.127	219.199.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.37mm	219.199.137	219.199.137S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.19mm	219.999.119	219.999.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.27mm	219.999.127	219.999.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	219.129.090	219.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.00mm	219.129.100	219.129.100S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.18mm	219.129.118	219.129.118S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	219.129.119	219.129.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.23mm	219.129.123	219.129.123S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	219.129.127	219.129.127S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.37mm	219.129.137	219.129.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	6/12.5mm	1.27mm	n/a	219.989.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12.5mm	1.27mm	219.127.127	n/a

Overview references

Length	Width	Thickness
Performance Blades		1.00mm
90mm	12mm	262.129.100 n/a

Note: sterile version with S

Performance Fits to: Aesculap®



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12mm	1.00mm	262.129.100	n/a

Overview references

Length	Width	Thickness				
		0.90mm	1.00mm	1.19mm	1.27mm	1.37mm
Performance Blades						
100mm	19mm				269.190.127 269.190.127S	
90mm	25mm		269.259.100 n/a		269.259.127 n/a	
90mm	19mm				269.199.127 269.199.127S	n/a 269.199.137S
90mm	12.5mm	269.129.090 269.129.090S	269.129.100 269.129.100S	269.129.119 n/a	269.129.127 269.129.127S	
80mm	12.5mm				269.128.127 269.128.127S	
80mm	8mm				269.088.127 269.088.127S	
70mm	25mm		269.257.100 n/a			

Note: sterile version with S

Performance Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.27mm	269.190.127	269.190.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	269.259.100	n/a

Performance Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	269.259.127	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	269.199.127	269.199.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.37mm	n/a	269.199.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	269.129.090	269.129.090S

Performance Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.00mm	269.129.100	269.129.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	269.129.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	269.129.127	269.129.127S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	12.5mm	1.27mm	269.128.127	269.128.127S

Performance Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	8mm	1.27mm	269.088.127	269.088.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	25mm	1.00mm	269.257.100	n/a

Overview references

Length	Width	Thickness				Description
		1.00mm	1.10mm	1.20mm	1.40mm	
Reciprocating Blades						
85mm	11.6mm			233.785.120 233.785.120S		Reciprocating aggressive
82mm	10mm	233.782.100 233.782.100S				Reciprocating aggressive
80mm	10mm				233.980.990 n/a	Reciprocating (file)
77mm	12mm		n/a 233.177.110S			Reciprocating
70mm	12mm		n/a 233.270.110S			Reciprocating (double-sided)
60mm	9mm	233.160.100 n/a				Reciprocating
34mm	7.2mm		233.334.110 233.334.110S			Reciprocating
cementless		n/a 233.018.275S				Reciprocating Keel Blade
cemented		n/a 233.018.290S				Reciprocating Keel Blade
cementless		n/a 233.000.001S				Reciprocating Keel Blade Set
cemented		n/a 233.000.002S				Reciprocating Keel Blade Set

Note: sterile version with S

Reciprocating Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	233.785.120	233.785.120S

Reciprocating Fits to: Acculan 3Ti




Length	Width	Thickness	REF non-sterile	REF sterile
82mm	10mm	1.00mm	233.782.100	233.782.100S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	10mm	1.4mm	233.980.990	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
77mm	12mm	1.10mm	n/a	233.177.110S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12mm	1.10mm	n/a	233.270.110S

Reciprocating Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
60mm	9mm	1.00mm	233.160.100	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
34mm	7.2mm	1.10mm	233.334.110	233.334.110S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	233.018.275S



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	233.018.290S

Keel Blade Set



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	219.129.090	219.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	233.785.120	233.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	233.018.275S

* Keel Blade Set only available sterile

Keel Blade Set cementless

REF non-sterile

REF sterile

n/a

233.000.001S

Keel Blade Set



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	219.129.090	219.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	233.785.120	233.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	233.018.290S

* Keel Blade Set only available sterile

Keel Blade Set cemented

REF non-sterile

REF sterile

n/a

233.000.002S

Overview references

Length	Width	Thickness
Trauma Blades		0.80mm
70mm	20mm	202.2070.080 n/a
48mm	14mm	202.1448.080 n/a

Note: sterile version with S

Trauma Fits to: Aesculap®



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	20mm	0.80mm	202.2070.080	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	14mm	0.80mm	202.1448.080	n/a

Overview references

Length	Width	Thickness		
		0.60mm	0.80mm	1.20mm
Trauma Blades				
70mm	27mm		209.2770.080 209.2770.080S	209.2770.120 209.2770.120S
70mm	20mm		209.2070.080 209.2070.080S	
70mm	18mm		209.1870.080 209.1870.080S	
50mm	10mm		209.1050.080 209.1050.080S	
48mm	27mm		209.2748.080 n/a	
48mm	14mm	209.1448.060 209.1448.060S		
35mm	10mm	209.1035.060 n/a	209.1035.080 209.1035.080S	
25mm	14mm	209.1425.060 209.1425.060S		
25mm	6mm	209.0625.060 209.0625.060S		

Note: sterile version with S

Trauma Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	0.80mm	209.2770.080	209.2770.080S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	1.20mm	209.2770.120	209.2770.120S

Trauma Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	20mm	0.80mm	209.2070.080	209.2070.080S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	18mm	0.80mm	209.1870.080	209.1870.080S



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	10mm	0.80mm	209.1050.080	209.1050.080S



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	27mm	0.60mm	209.2748.060	n/a

Trauma Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	14mm	0.60mm	209.1448.060	209.1448.060S



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	10mm	0.60mm	209.1035.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	10mm	0.80mm	209.1035.080	209.1035.080S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	14mm	0.60mm	209.1425.060	209.1425.060S

Trauma Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	209.1425.060	209.1425.060S

Overview references

Length	Width	Thickness	
		0.40mm	0.60mm
Small Bone Blades			
50mm	19mm		289.1950.060 289.1950.060S
50mm	10mm		289.1050.060 289.1050.060S
31mm	7mm	289.0731.040F 289.0731.040FS	
25mm	14mm		289.1425.060 289.1425.060S
25mm	10mm		289.1025.060 n/a
25mm	9mm	289.0925.040 289.0925.040S	
25mm	6mm		289.0625.060 289.0625.060S
21mm	9mm	289.0921.040F 289.0921.040FS	

Note: sterile version with S

Small Bone Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	19mm	0.60mm	289.1950.060	289.1950.060S



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	10mm	0.60mm	289.1050.060	289.1050.060S

Small Bone Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	7mm	0.40mm	289.0731.040F	289.0731.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	14mm	0.60mm	289.1425.060	289.1425.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	10mm	0.60mm	289.1025.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	289.0925.040	289.0925.040S

Small Bone Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	289.0625.060	289.0625.060S



Length	Width	Thickness	REF non-sterile	REF sterile
21mm	9mm	0.40mm	289.0921.040F	289.0921.040FS

Overview references

Length	Width	Thickness	
Small Bone Blades		0.40mm	0.60mm
41mm	14mm		295.0404.060 295.0404.060S
41mm	9mm		295.0403.060 n/a
25mm	9mm	295.0402.040 295.0402.040S	
25mm	6mm		295.0401.060 295.0401.060S

Note: sterile version with S

Small Bone Fits to: Elan 4, Acculan 4



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	14mm	0.60mm	295.0404.060	295.0404.060S



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	9mm	0.60mm	295.0403.060	n/a

Small Bone Fits to: Elan 4, Acculan 4



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	295.0402.040	295.0402.040S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	295.0401.060	295.0401.060S

About Gomina



Gomina is a leading manufacturer of medical products for bone surgery. We develop and produce our surgical instruments exclusively in Switzerland and market them worldwide. Hightech products from the heart of the Swiss mountains.

Gomina has pioneered the development of effective tooth geometry for oscillating saw blades using a combination of sharpened cutting and raker teeth. This unique combination of saw teeth results in outstanding cutting performance. Gomina's tooth geometry and arrangement minimises any lateral swerve during cutting, simplifying handling for the clinician.

Gomina's Keel Blade with its patented, innovative tooth geometry enables the accurate preparation of grooves without bone fragments.

100% SWISS MADE

The high requirements of saw technology commanded by surgeons form the basis for our innovation. Today our innovative medical products are used routinely in many knee, hip, foot and hand surgeries to the benefit of patients.

PRECISION AND SWISS-QUALITY

Since the founding of Gomina in 1974, we have made it our business to be precise.

Since then, Gomina has become a leader in the production of micro-mechanical precision components.

Precision that is used around the world in medical technology and sensor technology.


PRODUCTION

From the initial blank to the packaged finished unsterile product, Gomina carries out all processing steps in-house, except the sterilization, in compliance with the Swiss Quality Label.

R&D

Our engineering department converts visionary ideas into innovative products. Thanks to the latest technologies we are able to fulfill special customer requests. In close cooperation with the surgeons we are constantly expanding our product portfolio.

Certificates



In compliance with the Swiss Quality Label and the standards of the Medical Device Regulation established by the European Union, Gomina fulfills and ensures the strict requirements for safe and efficient medical products. All Gomina products are MDR-certified. Furthermore we fulfill the quality standards within the ISO 13485.

Our saw blades are certified for reusability, except for the Reciprocating Keel Blade. The valid certificates, instructions for use, processing and cleaning documents can be found online on www.gomina.ch

Compatibility

Gomina is a renowned manufacturer of high quality surgical saw blades for the use in modern orthopaedic, traumatological and reconstructive surgery. Our products are compatible with the power systems of various third party manufacturers established on the market. The compatibility of our replacement saw blades with the different power systems is indicated in this catalogue and on our packaging labels.

Our products are safe and effective to use in, but are not limited to, the power systems of Aesculap, Arthrex, Aygun, Conmed Linvatec Hall, DePuy Synthes, DeSoutter Medical, MatOrtho, MicroAire, Smith & Nephew, Bojin, Stryker and Zimmer Biomet. All trademarks mentioned in this document are registered trademarks of the respective companies. Gomina has no association with or is in any way affiliated with the mentioned companies. Data is correct at the time of publishing and its subject to change without notice.

In the past, uncertainty arose among some users, whether or not it is permitted to use Gomina saw blades (third party product) with their power systems, or if they are obliged to use the „original accessories“ provided by the manufacturer of the power system. This is due to ambiguous statements in the IFUs of some power systems or to misleading information given by the employees of individual power system manufacturers or distributors.

In this context, we would like to confirm that the use of Gomina's surgical saw blades in the power systems specified by us is definitely legally permitted and safe.

For any questions or additional information, please contact us.



Contact



Gomina AG
Raiftstrasse 4
3989 Niederwald, Switzerland
Phone +41 27 970 12 70
Fax +41 27 970 12 75
www.gomina.ch
sales@gomina.ch



GOMINA Deutschland GmbH
Ludwig-Vierling-Strasse 7
89233 Neu-Ulm, Germany
info@gomina.de

Version 2026-03
Last modification: 10.03.2026