

**Surgical Saw Blades for Orthopaedics,  
Trauma and Small Bones**  
PRODUCT CATALOGUE



REUSABLE  
SAW BLADES  
(except Keel Blade)



# Overview Gomina Saw Blades

## Saw blades for Orthopaedics



## Saw blades for Trauma



## Saw blades for Small Bones



## Reciprocating Saw blades



## Sternum Saw blades






























## Reciprocating Keel Blade



## Keel Blade Set



# Saw Blade Portfolio

						Page			
Fits to Aesculap®		Aesculap® Acculan 3TI GA673, GB660R, Acculan 4 GA331, GA336		Aesculap® Acculan 3TI GB435, GB440, GA635, GA636, GA641, GA674 Acculan 4 GA334		Aesculap® Elan 4 System GA673, GB660R, GB891R	10		
Fits to Arthrex®		Arthrex Power System 600™ AR-600SAG		Arthrex Power System 600™, DrillSaw Sport 400™ AR-600SAGMIS, AR-400DS		Arthrex Power System 600™, DrillSaw Mini 300™, DrillSaw Sport 400™ AR-300SAG, AR-400SAG, AR-600SAGS	38		
Fits to Aygun® Karınca 206		Aygun® Karınca 206 M206-300					49		
Fits to De Soutter®		De Soutter® Orthodrive, Multidrive DBK-700®, DPK-676®, IX-600®		De Soutter® Traumatdrive NQ700/707®, NU-470/2®, NZ-450®, NI-270®, NV-220®		De Soutter® Traumatdrive CU-471®, CZ-450®, CC-700/707®, DC-700®, CQ-700/707®,	50		
Fits to Linvatec® / Hall®		Linvatec®/ Hall® PowerPro®, Versipower®, Versipower Plus®, Series 3®, Series 4®		Linvatec®/ Hall® PowerPro®, Versipower®, Versipower Plus®, Series 3®, Series 4®		Linvatec®/ Hall® Versipower®, Versipower Plus®, Series 3®, Series 4®		Linvatec®/ Hall® MicroChose®, Micro100®, Micro E®	62
Fits to MicroAire®		MicroAire® Series 7000					91		
Fits to Stryker®		Stryker® EHD System 2000, System 4, System 5, System 6, System 7, System 8, System 9		Stryker® EHD System 2000, System 4, System 5, System 6, System 7, System 8, System 9		Stryker® TPS, TPX, F1	92		
Fits to Synthes® (AO)		Sodem™ Synthes® TRS BPL II CAD II		Synthes®		Synthes® Calibri, UNUM™		Synthes® E-Pen	123
Fits to Tuke Saw™		Tuke Saw™ System					157		
Fits to Zimmer®		Zimmer® X-Series Power System		Zimmer® Universal Power System		AO Universal Power System	159		
Fits to Bojin®		Bojin® System 1300					165		

# Teeth designs

## Premium teeth design

The Gomina Premium Line of saw blades is characterised by a combination of cross-grinded cutting teeth and straight-grinded raker teeth. This results in optimal cutting efficiency and high precision. The teeth are aligned on a straight line and are asymmetrical, which leads to a high degree of self-centring of the saw blade.



Premium teeth design

## Performance teeth design

The Gomina Performance Line of saw blades is characterised by asymmetrical straight-grinded cutting teeth, slightly less aggressive than the Premium Line.




Performance teeth design

## Reciprocating teeth designs

The Gomina reciprocating saw blades offers a variety of teeth designs. The cross-grinded cutting teeth (reciprocating aggressive) results in optimal cutting efficiency and high precision.

Reciprocating saw blades with the straight-grinded and symmetrical cutting teeth design are less aggressive. The set teeth allow an excellent sawing performance without jamming. Gomina double-sided reciprocating blades completes the portfolio.




## Trauma teeth design

The Gomina trauma saw blades offers the wide product range for the corresponding application. The teeth set allow an excellent sawing performance without jamming.




## Small Bone teeth design

The Gomina small bone saw blade is characterised by asymmetrical straight-grinded cutting teeth. These saw blades are mostly used without cutting blocks due the ease for free-hand sawing.



# General Saw Blades information



NOTE: All markings on the products in this catalogue are exemplary, positions may differ from the actual marking. The UDI is to be understood as an example. Every article has its individual code.

# Packaging

## NON-STERILE PACKAGING

Non-sterile saw blades are packaged in single pouches and carry the product label with detailed information (label printed black & white). They don't come in boxes and can be **ordered single** for any desired quantity. Our unsterile saw blades are all reusable.



## STERILE PACKAGING

Sterile saw blades are packaged in double sterile pouches which feature an additional protection layer for the cutting edge. Minimal ordering quantity is **one box containing 5 saw blades**. The box also contains one patient label per saw blade, placed between the first and second pouch of each blade. Except for the Reciprocating Keel Blade, all our sterile saw blades are reusable.



# Contents

<b>Fits to Aesculap®</b>	<b>10</b>
Premium Fits to: Acculan 3Ti	10
Performance Fits to: Aesculap®	17
Performance Fits to: Acculan 3Ti	18
Reciprocating Fits to: Acculan 3Ti	22
Keel Blade Set	25
Trauma Fits to: Aesculap®	27
Trauma Fits to: Acculan 3Ti	28
Small Bone Fits to: Acculan 3Ti	32
Small Bone Fits to: Elan 4, Acculan 4	35
<b>Fits to Arthrex®</b>	<b>37</b>
Premium Fits to: Power System 600	37
Premium Fits to: AR-600-SAGMIS, AR-400DS	39
Performance Fits to: AR-600-SAGMIS, AR-400DS	41
Small Bone Fits to: AR-600-SAGMIS, AR-400DS	42
Small Bone	44
<b>Fits to Aygun® Karınca 206</b>	<b>48</b>
Premium	48
<b>Fits to De Soutter®</b>	<b>49</b>
Premium Fits to: Multidrive™, Orthodrive®	49
Performance Fits to: Multidrive™, Orthodrive®	52
Reciprocating Fits to: Multidrive™	54
Reciprocating Fits to: Orthodrive® DBC-700	55
Keel Blade Set	56
Trauma Fits to: Multidrive™, Orthodrive®	58
Small Bone Fits to: Traumadrive™	59
<b>Fits to Linvatec® / Hall®</b>	<b>61</b>
Premium Fits to: Power Pro®, Versipower®Plus	61
Premium Fits to: Series 3 & 4, Versipower	68
Performance Fits to: Power Pro®, Versipower®Plus	69
Performance Fits to: Series 3&4, Versipower	74
Reciprocating	75
Keel Blade Set	78

Trauma Fits to: Power Pro®, Versipower®Plus	80
Small Bone (new hub)	82
Small Bone (old hub)	86
Small Bone	88
<b>Fits to MicroAire®</b>	<b>90</b>
Performance	90
<b>Fits to Stryker®</b>	<b>91</b>
Premium	91
Performance	101
Reciprocating	107
Keel Blade Set	111
Trauma	113
Small Bone	116
<b>Fits to Synthes® (AO)</b>	<b>122</b>
Premium	122
Performance	130
Reciprocating	135
Keel Blade Set	139
Trauma	141
Small Bones Fits to: Colibri, UNIUM™	146
Small Bones Fits to: Electric Pen Drive	152
<b>Fits to Tuke Saw™</b>	<b>156</b>
Reciprocating	156
<b>Fits to Zimmer®</b>	<b>158</b>
Performance Fits to: X-Series	158
Reciprocating	163
<b>Fits to Bojin®</b>	<b>164</b>
Premium	164

# Fits to Aesculap®

## Overview references

Length	Width	Thickness						
		0.90mm	1.00mm	1.18mm	1.19mm	1.23mm	1.27mm	1.37mm
<b>Premium Blades</b>								
100mm	19mm				219.190.119 219.190.119S			
90mm	25mm	219.259.090 219.259.090S	219.259.100 219.259.100S	219.259.118 219.259.118S	219.259.119 219.259.119S		219.259.127 219.259.127S	219.259.137 219.259.137S
90mm	19/22mm				219.969.119 219.969.119S		219.969.127 219.969.127S	
90mm	19mm	219.199.090 219.199.090S			219.199.119 219.199.119S	219.199.123 n/a	219.199.127 219.199.127S	219.199.137 219.199.137S
90mm	12.5/19mm				219.999.119 219.999.119S		219.999.127 219.999.127S	
90mm	12.5mm	219.129.090 219.129.090S	219.129.100 219.129.100S	219.129.118 219.129.118S	219.129.119 219.129.119S	219.129.123 219.129.123S	219.129.127 219.129.127S	219.129.137 219.129.137S
90mm	6/12.5mm						n/a 219.989.127S	
70mm	12.5mm						219.127.127 n/a	

Note: sterile version with S

## Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.19mm	219.190.119	219.190.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	0.90mm	219.259.090	219.259.090S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	219.259.100	219.259.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.18mm	219.259.118	219.259.118S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.19mm	219.259.119	219.259.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	219.259.127	219.259.127S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.37mm	219.259.137	219.259.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19/22mm	1.19mm	219.969.119	219.969.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19/22mm	1.27mm	219.969.127	219.969.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	0.90mm	219.199.090	219.199.090S

Premium Fits to Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	219.199.119	219.199.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.23mm	219.199.123	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	219.199.127	219.199.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.37mm	219.199.137	219.199.137S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.19mm	219.999.119	219.999.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.27mm	219.999.127	219.999.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	219.129.090	219.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.00mm	219.129.100	219.129.100S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.18mm	219.129.118	219.129.118S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	219.129.119	219.129.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.23mm	219.129.123	219.129.123S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	219.129.127	219.129.127S

Premium Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.37mm	219.129.137	219.129.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	6/12.5mm	1.27mm	n/a	219.989.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12.5mm	1.27mm	219.127.127	n/a

Overview references

Length	Width	Thickness
<b>Performance Blades</b>		<b>1.00mm</b>
90mm	12mm	262.129.100 n/a

Note: sterile version with S

Performance Fits to: Aesculap®



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12mm	1.00mm	262.129.100	n/a

Overview references

Length	Width	Thickness				
		0.90mm	1.00mm	1.19mm	1.27mm	1.37mm
<b>Performance Blades</b>						
100mm	19mm				269.190.127 269.190.127S	
90mm	25mm		269.259.100 n/a		269.259.127 n/a	
90mm	19mm				269.199.127 269.199.127S	n/a 269.199.137S
90mm	12.5mm	269.129.090 269.129.090S	269.129.100 269.129.100S	269.129.119 n/a	269.129.127 269.129.127S	
80mm	12.5mm				269.128.127 269.128.127S	
80mm	8mm				269.088.127 269.088.127S	
70mm	25mm		269.257.100 n/a			

Note: sterile version with S

Performance Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.27mm	269.190.127	269.190.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	269.259.100	n/a

Performance Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	269.259.127	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	269.199.127	269.199.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.37mm	n/a	269.199.137S




Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	269.129.090	269.129.090S

Performance Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.00mm	269.129.100	269.129.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	269.129.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	269.129.127	269.129.127S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	12.5mm	1.27mm	269.128.127	269.128.127S

Performance Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	8mm	1.27mm	269.088.127	269.088.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	25mm	1.00mm	269.257.100	n/a

Overview references

Length	Width	Thickness				Description
		1.00mm	1.10mm	1.20mm	1.40mm	
<b>Reciprocating Blades</b>						
85mm	11.6mm			233.785.120 233.785.120S		Reciprocating aggressive
82mm	10mm	233.782.100 233.782.100S				Reciprocating aggressive
80mm	10mm				233.980.990 n/a	Reciprocating (file)
77mm	12mm		n/a 233.177.110S			Reciprocating
70mm	12mm		n/a 233.270.110S			Reciprocating (double-sided)
60mm	9mm	233.160.100 n/a				Reciprocating
34mm	7.2mm		233.334.110 233.334.110S			Reciprocating
cementless		n/a 233.018.275S				Reciprocating Keel Blade
cemented		n/a 233.018.290S				Reciprocating Keel Blade
cementless		n/a 233.000.001S				Reciprocating Keel Blade Set
cemented		n/a 233.000.002S				Reciprocating Keel Blade Set

Note: sterile version with S

Reciprocating Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	233.785.120	233.785.120S

Reciprocating Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
82mm	10mm	1.00mm	233.782.100	233.782.100S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	10mm	1.4mm	233.980.990	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
77mm	12mm	1.10mm	n/a	233.177.110S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12mm	1.10mm	n/a	233.270.110S

Reciprocating Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
60mm	9mm	1.00mm	233.160.100	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
34mm	7.2mm	1.10mm	233.334.110	233.334.110S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	233.018.275S



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	233.018.290S

**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	219.129.090	219.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	233.785.120	233.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	233.018.275S

\* Keel Blade Set only available sterile

**Keel Blade Set cementless**

**REF non-sterile**

**REF sterile**

n/a

233.000.001S

**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	219.129.090	219.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	233.785.120	233.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	233.018.290S

\* Keel Blade Set only available sterile

**Keel Blade Set cemented**

**REF non-sterile**

**REF sterile**

n/a

233.000.002S

Overview references

Length	Width	Thickness
<b>Trauma Blades</b>		<b>0.80mm</b>
70mm	20mm	202.2070.080 n/a
48mm	14mm	202.1448.080 n/a

Note: sterile version with S

Trauma Fits to: Aesculap®



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	20mm	0.80mm	202.2070.080	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	14mm	0.80mm	202.1448.080	n/a

Overview references

Length	Width	Thickness		
		0.60mm	0.80mm	1.20mm
<b>Trauma Blades</b>				
70mm	27mm		209.2770.080 209.2770.080S	209.2770.120 209.2770.120S
70mm	20mm		209.2070.080 209.2070.080S	
70mm	18mm		209.1870.080 209.1870.080S	
50mm	10mm		209.1050.080 209.1050.080S	
48mm	27mm		209.2748.080 n/a	
48mm	14mm	209.1448.060 209.1448.060S		
35mm	10mm	209.1035.060 n/a	209.1035.080 209.1035.080S	
25mm	14mm	209.1425.060 209.1425.060S		
25mm	6mm	209.0625.060 209.0625.060S		

Note: sterile version with S

Trauma Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	0.80mm	209.2770.080	209.2770.080S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	1.20mm	209.2770.120	209.2770.120S

Trauma Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	20mm	0.80mm	209.2070.080	209.2070.080S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	18mm	0.80mm	209.1870.080	209.1870.080S




Length	Width	Thickness	REF non-sterile	REF sterile
50mm	10mm	0.80mm	209.1050.080	209.1050.080S



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	27mm	0.60mm	209.2748.060	n/a

Trauma Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	14mm	0.60mm	209.1448.060	209.1448.060S



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	10mm	0.60mm	209.1035.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	10mm	0.80mm	209.1035.080	209.1035.080S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	14mm	0.60mm	209.1425.060	209.1425.060S

Trauma Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	209.1425.060	209.1425.060S

Overview references

Length	Width	Thickness	
		0.40mm	0.60mm
<b>Small Bone Blades</b>			
50mm	19mm		289.1950.060 289.1950.060S
50mm	10mm		289.1050.060 289.1050.060S
31mm	7mm	289.0731.040F 289.0731.040FS	
25mm	14mm		289.1425.060 289.1425.060S
25mm	10mm		289.1025.060 n/a
25mm	9mm	289.0925.040 289.0925.040S	
25mm	6mm		289.0625.060 289.0625.060S
21mm	9mm	289.0921.040F 289.0921.040FS	

Note: sterile version with S

Small Bone Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	19mm	0.60mm	289.1950.060	289.1950.060S



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	10mm	0.60mm	289.1050.060	289.1050.060S

Small Bone Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	7mm	0.40mm	289.0731.040F	289.0731.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	14mm	0.60mm	289.1425.060	289.1425.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	10mm	0.60mm	289.1025.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	289.0925.040	289.0925.040S

Small Bone Fits to: Acculan 3Ti



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	289.0625.060	289.0625.060S



Length	Width	Thickness	REF non-sterile	REF sterile
21mm	9mm	0.40mm	289.0921.040F	289.0921.040FS

Overview references

Length	Width	Thickness	
<b>Small Bone Blades</b>		<b>0.40mm</b>	<b>0.60mm</b>
41mm	14mm		295.0404.060 295.0404.060S
41mm	9mm		295.0403.060 n/a
25mm	9mm	295.0402.040 295.0402.040S	
25mm	6mm		295.0401.060 295.0401.060S

Note: sterile version with S

Small Bone Fits to: Elan 4, Acculan 4



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	14mm	0.60mm	295.0404.060	295.0404.060S



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	9mm	0.60mm	295.0403.060	n/a

Small Bone Fits to: Elan 4, Acculan 4



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	295.0402.040	295.0402.040S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	295.0401.060	295.0401.060S

# Fits to Arthrex®

## Overview references

Length	Width	Thickness	
<b>Premium Blades</b>		<b>0.90mm</b>	<b>1.27mm</b>
90mm	19mm		221.191.127 221.191.127S
90mm	19mm		221.199.127 221.199.127S
90mm	12.5mm	221.129.090 n/a	221.129.127 221.129.127S

Note: sterile version with S

### Premium Fits to: Power System 600



Length	Width	Thickness	REF non-sterile	REF sterile
110mm	19mm	1.27mm	221.191.127	221.191.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	221.199.127	221.199.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	221.129.090	n/a

Premium Fits to: Power System 600



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	221.129.127	221.129.127S

# Fits to Arthrex®

## Overview references

Length	Width	Thickness		
Premium Blades		0.90mm	1.19mm	1.27mm
90mm	25mm	294.0205.090 n/a		
90mm	19mm	294.0202.090 n/a		
90mm	19mm		294.0203.119 n/a	294.0203.127 294.0203.127S
90mm	12.5mm	294.0201.090 n/a		

Note: sterile version with S

## Premium Fits to: AR-600-SAGMIS, AR-400DS



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	0.90mm	294.0205.090	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	0.90mm	294.0202.090	n/a


Premium Fits to: AR-600-SAGMIS, AR-400DS



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	294.0203.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	294.0203.127	294.0203.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	294.0201.090	n/a

# Fits to Arthrex®

## Overview references

Length	Width	Thickness
<b>Performance Blades</b>		<b>1.27mm</b>
80mm	12.5mm	n/a 294.0206.127S
80mm	8mm	n/a 294.0207.127S

Note: sterile version with S

## Performance Fits to: AR-600-SAGMIS, AR-400DS



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	12.5mm	1.27mm	n/a	294.0206.127S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	8mm	1.27mm	n/a	294.0207.127S

# Fits to Arthrex®

## Overview references

Length	Width	Thickness	
<b>Small Bone Blades</b>		<b>0.40mm</b>	<b>0.60mm</b>
70mm	19mm		295.0303.060 n/a
50mm	19mm		295.0302.060 n/a
45mm	9mm		295.0304.060 n/a
31mm	9mm	295.0305.040 n/a	
25mm	9mm	295.0306.040 n/a	
25mm	9mm	295.0307.040 n/a	

Note: sterile version with S

## Small Bone Fits to: AR-600-SAGMIS, AR-400DS



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	0.60mm	295.0303.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	19mm	0.60mm	295.0302.060	n/a

**Small Bone Fits to: AR-600-SAGMIS, AR-400DS**



Length	Width	Thickness	REF non-sterile	REF sterile
45mm	9mm	0.60mm	295.0304.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	9mm	0.40mm	295.0305.040	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	295.0306.040	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.40mm	295.0307.040	n/a

Overview references

Length	Width	Thickness	
<b>Small Bone Blades</b>			
		<b>0.40mm</b>	<b>0.60mm</b>
70mm	19mm		287.1970.060 n/a
50mm	27mm		287.2750.060 n/a
50mm	19mm		287.1950.060 n/a
41mm	14mm		287.1441.060 287.1441.060S
41mm	9mm		287.0941.060 287.0941.060S
41mm	9mm		287.0941.060F 287.0941.060FS
35mm	9mm		287.0935.060 n/a
25mm	9mm	287.0925.040 287.0925.040S	287.0925.060F 287.0925.060FS
25mm	9mm	287.0925.040F 287.0925.040FS	
25mm	6mm		287.0625.060 287.0625.060S
18mm	5mm	287.0518.040F 287.0518.040FS	

Note: sterile version with S

Small Bone



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	0.60mm	287.1970.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	27mm	0.60mm	287.2750.060	n/a

Small Bone



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	19mm	0.60mm	287.1950.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	14mm	0.60mm	287.1441.060	287.1441.060S



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	9mm	0.60mm	287.0941.060	287.0941.060S



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	9mm	0.60mm	287.0941.060F	287.0941.060FS

Small Bone



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	9mm	0.60mm	287.0935.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	287.0925.040	287.0925.040S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	287.0925.040F	287.0925.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.60mm	287.0925.060F	287.0925.060FS

**Small Bone**



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	287.0625.060	287.0625.060S



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	5mm	0.40mm	287.0518.040F	287.0518.040FS

# Fits to Aygun® Karınca 206

## Overview references

Length	Width	Thickness
<b>Premium Blades</b>		<b>1.27mm</b>
100mm	19mm	294.0301.127 n/a
100mm	12/19mm	294.0302.127 n/a

Note: sterile version with S

## Premium



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.27mm	294.0301.127	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	12/19mm	1.27mm	294.0302.127	n/a

# Fits to De Soutter®

## Overview references

Length	Width	Thickness				
		0.90mm	1.19mm	1.23mm	1.27mm	1.37mm
<b>Premium Blades</b>						
90mm	25mm			n/a 217.259.123S	217.259.127 217.259.127S	
90mm	19mm				217.199.127 217.199.127S	217.199.137 n/a
90mm	12.5mm	217.129.090 217.129.090S	217.129.119 217.129.119S	217.129.123 217.129.123S		

Note: sterile version with S

## Premium Fits to: Multidrive™, Orthodrive®



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.23mm	n/a	217.259.123S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27	217.259.127	217.259.127S

Premium Fits to: Multidrive™, Orthodrive®



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	217.199.127	217.199.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.37mm	217.199.137	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	217.129.090	217.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	217.129.119	217.129.119S

Premium Fits to: Multidrive™, Orthodrive®



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.23mm	217.129.123	217.129.123S

Overview references

Length	Width	Thickness				
		0.90mm	1.00mm	1.19mm	1.27mm	1.37mm
<b>Performance Blades</b>						
90mm	25mm		267.259.100 267.259.100S			267.259.137 267.259.137S
90mm	19mm	267.199.090 n/a		267.199.119 267.199.119S		
80mm	12.5mm				267.128.127 267.128.127S	
80mm	8mm				267.088.127 267.088.127S	

Note: sterile version with S

Performance Fits to: Multidrive™, Orthodrive®



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	267.259.100	267.259.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.37mm	267.259.137	267.259.137S

Performance Fits to: Multidrive™, Orthodrive®



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	0.90mm	267.199.090	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	267.199.119	267.199.119S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	12.5mm	1.27mm	267.128.127	267.128.127S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	8mm	1.27mm	267.088.127	267.088.127S

Overview references

Length	Width	Thickness		Description
<b>Reciprocating Blades</b>				
		1.10mm	1.20mm	
85mm	11.6mm		231.785.120 231.785.120S	Reciprocating aggressive (Multidrive)
85mm	11.6mm		296.0504.120 296.0504.120S	Reciprocating (Orthodrive DBC-700)
77mm	12mm	296.0505.110 296.0505.110S		Reciprocating (Orthodrive DBC-700)
cementless		n/a 231.018.275S		Reciprocating Keel Blade
cemented		n/a 231.018.290S		Reciprocating Keel Blade
cementless		n/a 231.000.003S		Reciprocating Keel Blade Set
cemented		n/a 231.000.004S		Reciprocating Keel Blade Set

Note: sterile version with S

Reciprocating Fits to: Multidrive™



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	231.785.120	231.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	231.018.275S

**Reciprocating Fits to: Multidrive™**



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	231.018.290S

**Reciprocating Fits to: Orthodrive® DBC-700**



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	296.0504.120	296.0504.120S



Length	Width	Thickness	REF non-sterile	REF sterile
77mm	12mm	1.10mm	296.0505.110	296.0505.110S

**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	217.129.090	217.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	231.785.120	231.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	231.018.275S

\* Keel Blade Set only available sterile

**Keel Blade Set cementless**

**REF non-sterile**

**REF sterile**

n/a

231.000.003S

**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	217.129.090	217.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	231.785.120	231.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	231.018.290S

\* Keel Blade Set only available sterile

**Keel Blade Set cemented**

**REF non-sterile**

**REF sterile**

n/a

231.000.004S

Overview references

Length	Width	Thickness	
<b>Trauma Blades</b>		<b>0.60mm</b>	<b>1.00mm</b>
64mm	34mm	207.3464.100 207.3464.100S	

Note: sterile version with S

Trauma Fits to: Multidrive™, Orthodrive®



Length	Width	Thickness	REF non-sterile	REF sterile
64mm	34mm	1.00mm	207.3464.100	207.3464.100S

Overview references

Length	Width	Thickness
<b>Small Bone Blades</b>		<b>0.60mm</b>
45mm	15mm	298.1545.060F 298.1545.060FS
25mm	10mm	298.1025.060 298.1025.060S
25mm	9mm	298.0925.060 298.0925.060S
22mm	6mm	298.0622.060 298.0622.060S
21mm	7mm	298.0721.060 298.0721.060S
18mm	9mm	298.0918.060 298.0918.060S

Note: sterile version with S

Small Bone Fits to: Traumadrive™



Length	Width	Thickness	REF non-sterile	REF sterile
45mm	15mm	0.60mm	298.1545.060F	298.1545.060FS



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	10mm	0.60mm	298.1025.060	298.1025.060S

Small Bone Fits to: Traumadrive™



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.60mm	298.0925.060	298.0925.060S



Length	Width	Thickness	REF non-sterile	REF sterile
22mm	6mm	0.60mm	298.0622.060	298.0622.060S



Length	Width	Thickness	REF non-sterile	REF sterile
21mm	7mm	0.60mm	298.0721.060	298.0721.060S



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	9mm	0.60mm	298.0918.060	298.0918.060S

# Fits to Linvatec® / Hall®

## Overview references

Length	Width	Thickness						
		0.90mm	1.00mm	1.19mm	1.23mm	1.27mm	1.37mm	1.47mm
<b>Premium Blades</b>								
100mm	19mm					213.190.127 213.190.127S		
90mm	25mm	213.259.090 213.259.090S		213.259.119 213.259.119S		213.259.127 213.259.127S	213.259.137 213.259.137S	213.259.147 n/a
90mm	19mm	213.199.090 213.199.090S		213.199.119 213.199.119S		213.199.127 213.199.127S	213.199.137 213.199.137S	
90mm	12.5/19mm			213.999.119 213.999.119S		213.999.127 213.999.127S		
90mm	12.5mm	213.129.090 213.129.090S		213.129.119 213.129.119S	213.129.123 213.129.123S	213.129.127 213.129.127S	213.129.137 213.129.137S	213.129.147 213.129.147S
90mm	6/12.5mm			213.989.119 n/a				
70mm	8/12.5mm			n/a 213.977.119S				
70mm	19mm	213.197.090 213.197.090S		213.197.119 213.197.119S		213.197.127 213.197.127S		
70mm	12.5mm		213.127.100 213.127.100S					

Note: sterile version with S

## Premium Fits to: Power Pro®, Versipower® Plus



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.27mm	213.190.127	213.190.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	0.90mm	213.259.090	213.259.090S

Premium Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.19mm	213.259.119	213.259.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	213.259.127	213.259.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.37mm	213.259.137	213.259.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.47mm	219.259.147	n/a

Premium Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	0.90mm	213.199.090	213.199.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	213.199.119	213.199.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	213.199.127	213.199.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.37mm	213.199.137	213.199.137S

Premium Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.19mm	213.999.119	213.999.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.27mm	213.999.127	213.999.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	213.129.090	213.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	213.129.119	213.129.119S

Premium Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.23mm	213.129.123	213.129.123S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	213.129.127	213.129.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.37mm	213.129.137	213.129.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.47mm	213.129.147	213.129.147S

Premium Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	6/12.5mm	1.19mm	213.989.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	8/12.5mm	1.19mm	n/a	213.977.119S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	0.90mm	213.197.090	213.197.090S




Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	1.19mm	213.197.119	213.197.119S

Premium Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	1.27mm	213.197.127	213.197.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12.5mm	1.00mm	213.127.100	213.127.100S

Overview references

Length	Width	Thickness
<b>Premium Blades</b>		<b>1.23mm</b>
90mm	19mm	293.199.123 n/a
90mm	12.5mm	293.129.123 n/a

Note: sterile version with S

Premium Fits to: Series 3 & 4, Versipower



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.23mm	293.199.123	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.23mm	293.129.123	n/a

Overview references

Length	Width	Thickness			
Performance Blades		0.90mm	1.00mm	1.19mm	1.27mm
100mm	19mm				263.191.127 263.191.127S
90mm	25mm		263.259.100 263.259.100S		263.259.127 263.259.127S
90mm	19mm			263.199.119 263.199.119S	263.199.127 263.199.127S
90mm	12.5/19mm				263.999.127 263.999.127S
90mm	12.5mm	263.129.090 263.129.090S	263.129.100 263.129.100S	263.129.119 263.129.119S	263.129.127 263.129.127S
80mm	12.5mm				n/a 263.128.127S
80mm	8mm				n/a 263.088.127S
70mm	19mm			263.197.119 263.197.119S	
70mm	12.5mm	263.127.090 263.127.090S		263.127.119 n/a	263.127.127 263.127.127S
60mm	10mm				263.106.127 263.106.127S

Note: sterile version with S

Performance Fits to: Power Pro®, Versipower® Plus



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.27mm	263.191.127	263.191.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	263.259.100	263.259.100S

Performance Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	263.259.127	263.259.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	263.199.119	263.199.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	263.199.127	263.199.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.27mm	263.999.127	263.999.127S

Performance Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	263.129.090	263.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.00mm	263.129.100	263.129.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	263.129.119	263.129.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	263.129.127	263.129.127S

Performance Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	12.5mm	1.27mm	n/a	263.128.127S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	8mm	1.27mm	n/a	263.088.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	1.19mm	263.197.119	263.197.119S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12.5mm	0.90mm	263.127.090	263.127.090S

Performance Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12.5mm	1.19mm	263.127.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12.5mm	1.27mm	263.127.127	263.127.127S



Length	Width	Thickness	REF non-sterile	REF sterile
60mm	10mm	1.27mm	263.106.127	263.106.127S

**Overview references**

Length	Width	Thickness
<b>Performance Blades</b>		<b>1.27mm</b>
70mm	19mm	294.0101.127 294.0101.127S

Note: sterile version with S

**Performance Fits to: Series 3&4, Versipower**



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	1.27mm	294.0101.127	294.0101.127S

### Overview references

Length	Width	Thickness			Description
		1.00mm	1.10mm	1.20mm	
<b>Reciprocating Blades</b>					
85mm	11.6mm			231.785.120 231.785.120S	Reciprocating aggressive
77mm	12mm		231.177.110 231.177.110S		Reciprocating
70mm	12mm		231.270.110 231.270.110S		Reciprocating
60mm	9mm	231.160.100 231.160.100S			Reciprocating
34mm	9.5mm	231.334.100 231.334.100S			Reciprocating (sternum)
cementless		n/a 231.018.275S			Reciprocating Keel Blade
cemented		n/a 231.018.290S			Reciprocating Keel Blade
cementless		n/a 231.000.001S			Reciprocating Keel Blade Set
cemented		n/a 231.000.002S			Reciprocating Keel Blade Set

Note: sterile version with S

### Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	231.785.120	231.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
77mm	12mm	1.10mm	231.177.110	231.177.110S

**Reciprocating**



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12mm	1.10mm	231.270.110	231.270.110S



Length	Width	Thickness	REF non-sterile	REF sterile
60mm	9mm	1.00mm	231.160.100	231.160.100S



Length	Width	Thickness	REF non-sterile	REF sterile
34mm	9.5mm	1.00mm	231.334.100	231.334.100S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	231.018.275S

**Reciprocating**



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	231.018.290S

**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	213.129.090	213.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	231.785.120	231.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	231.018.275S

\* Keel Blade Set only available sterile

**Keel Blade Set cementless**

**REF non-sterile**

**REF sterile**

n/a

231.000.001S

**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	213.129.090	213.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	231.785.120	231.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	231.018.290S

\* Keel Blade Set only available sterile

**Keel Blade Set cemented**

**REF non-sterile**

**REF sterile**

n/a


231.000.002S

Overview references

Length	Width	Thickness			
		0.60mm	0.80mm	0.90mm	1.00mm
<b>Trauma Blades</b>					
70mm	27mm		203.2770.080 203.2770.080S		
70mm	20mm		203.2070.080 203.2070.080S		
64mm	34mm				203.3464.100 203.3464.100S
63mm	49mm			203.4963.090 n/a	
35mm	10mm	203.1035.060 203.1035.060S			
25mm	10mm	203.1025.060 203.1025.060S			

Note: sterile version with S

Trauma Fits to: Power Pro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	0.80mm	203.2770.080	203.2770.080S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	20mm	0.80mm	203.2070.080	203.2070.080S

Trauma Fits to: PowerPro®, Versipower®Plus



Length	Width	Thickness	REF non-sterile	REF sterile
64mm	34mm	1.00mm	203.3464.100	203.3464.100S



Length	Width	Thickness	REF non-sterile	REF sterile
63mm	49mm	0.90mm	203.4963.090	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	10mm	0.60mm	203.1035.060	203.1035.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	10mm	0.60mm	203.1025.060	203.1025.060S

Overview references

Length	Width	Thickness	
Small Bone Blades		0.40mm	0.60mm
41mm	19mm		283.1941.060 283.1941.060S
41mm	14mm		283.1441.060 283.1441.060S
41mm	9mm		283.0941.060 283.0941.060S
35mm	9mm		283.0935.060 283.0935.060S
25mm	14mm		283.1425.060 283.1425.060S
25mm	9mm	283.0925.040 283.0925.040S	283.0925.060 283.0925.060S
25mm	9mm	283.0925.040F 283.0925.040FS	
25mm	6mm	283.0625.040 283.0625.040S	283.0625.060 283.0625.060S
18mm	6mm	283.0618.040 283.0618.040S	
16.4mm	14mm		283.1001.060 n/a
12mm	9mm		283.0999.060 283.0999.060S

Note: sterile version with S

Small Bone (new hub)



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	19mm	0.60mm	283.1941.060	283.1941.060S



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	14mm	0.60mm	283.1441.060	283.1441.060S

**Small Bone (new Hub)**



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	9mm	0.60mm	283.0941.060	283.0941.060S



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	9mm	0.60mm	283.0935.060	283.0935.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	14mm	0.60	283.1425.060	283.1425.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	283.0925.040	283.0925.040S

**Small Bone (new Hub)**



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.60mm	283.0925.060	283.0925.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	283.0925.040F	283.0925.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.40mm	283.0625.040	283.0625.040S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	283.0625.060	283.0625.060S

**Small Bone (new Hub)**



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	6mm	0.40mm	283.0618.040	283.0618.040S



Length	Width	Thickness	REF non-sterile	REF sterile
16.4mm	14mm	0.60mm	283.1001.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
12mm	9mm	0.60mm	283.0999.060	283.0999.060S

### Overview references

Length	Width	Thickness
<b>Small Bone Blades</b>		<b>0.40mm</b>
35mm	11mm	284.1135.040 284.1135.040S
25mm	5mm	284.0525.040 284.0525.040S
18mm	9mm	284.0918.040 284.0918.040S
18mm	5mm	284.0518.040 284.0518.040S
14mm	5mm	284.0514.040 284.0514.040S

Note: sterile version with S

### Small Bone (old hub)



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	19mm	0.40mm	284.1135.040	284.1135.040S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	5mm	0.40mm	284.0525.040	284.0525.040S

**Small Bone (old Hub)**



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	9mm	0.40mm	284.0918.040	284.0918.040S



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	5mm	0.40mm	284.0518.040	284.0518.040S



Length	Width	Thickness	REF non-sterile	REF sterile
14mm	5mm	0.40mm	284.0514.040	284.0514.040S

### Overview references

Length	Width	Thickness
<b>Small Bone Blades</b>		<b>0.60mm</b>
41mm	14mm	295.0505.060 n/a
40mm	9mm	295.0502.060 295.0502.060S
35mm	9mm	295.0504.060 n/a
25mm	9mm	295.0503.060 295.0503.060S

Note: sterile version with S

### Small Bone




Length	Width	Thickness	REF non-sterile	REF sterile
41mm	14mm	0.60mm	295.0505.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
40mm	9mm	0.60mm	295.0502.060	295.0502.060S

**Small Bone**



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	9mm	0.60mm	295.0504.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.60mm	295.0503.060	295.0503.060S

# Fits to MicroAire®

## Overview references

Length	Width	Thickness
<b>Performance Blades</b>		<b>1.27mm</b>
90mm	25mm	270.259.127 270.259.127S

Note: sterile version with S

## Performance



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	270.259.127	270.259.127S

# Fits to Stryker®

## Overview references

Length	Width	Thickness							
Premium Blades		0.90mm	1.00mm	1.19mm	1.23mm	1.25mm	1.27mm	1.37mm	1.47mm
110mm	25mm						218.251.127 218.251.127S		
110mm	19mm						218.191.127 218.191.127S		
100mm	25mm			218.250.119 n/a			218.250.127 n/a		
100mm	19mm			218.190.119 218.190.119S			218.190.127 218.190.127S		
100mm	12.5/19mm						218.990.127 218.990.127S		
90mm	25mm	218.259.090 218.259.090S	218.259.100 218.259.100S	218.259.119 218.259.119S	218.259.123 218.259.123S	218.259.125 218.259.125S	218.259.127 218.259.127S	218.259.137 218.259.137S	218.259.147 218.259.147S
90mm	22mm							n/a 218.229.137S	
90mm	19mm	218.199.090 218.199.090S	218.199.100 218.199.100S	218.199.119 218.199.119S	218.199.123 218.199.123S		218.199.127 218.199.127S	218.199.137 218.199.137S	
90mm	12.5/19mm			218.999.119 218.999.119S			218.999.127 218.999.127S		
90mm	12.5mm	218.129.090 218.129.090S	218.129.100 n/a	218.129.119 218.129.119S	218.129.123 218.129.123S		218.129.127 218.129.127S	218.129.137 218.129.137S	
70mm	25mm	218.257.090 218.257.090S							
70mm	19mm						218.197.127 218.197.127S		218.197.147 218.197.147S
70mm	6/12.5mm						218.987.127 218.987.127S		218.987.147 218.987.147S
60mm	12.5mm	218.126.090 218.126.090S							

Note: sterile version with S

## Premium



Length	Width	Thickness	REF non-sterile	REF sterile
110mm	25mm	1.27mm	218.251.127	218.251.127S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
110mm	19mm	1.27mm	218.191.127	218.191.127S



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	25mm	1.19mm	218.250.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	25mm	1.27mm	218.250.127	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.19mm	218.190.119	218.190.119S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.27mm	218.190.127	218.190.127S



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	12.5/19mm	1.27mm	218.990.127	218.990.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	0.90mm	218.259.090	218.259.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	218.259.100	218.259.100S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.19mm	218.259.119	218.259.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.23mm	218.259.123	218.259.123S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.25mm	218.259.125	218.259.125S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	218.259.127	218.259.127S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.37mm	218.259.137	218.259.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.47mm	218.259.147	218.259.147S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	22mm	1.37mm	n/a	218.229.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	0.90mm	218.199.090	218.199.090S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.00mm	218.199.100	218.199.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	218.199.119	218.199.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.23mm	218.199.123	218.199.123S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	218.199.127	218.199.127S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.37mm	218.199.137	218.199.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.19mm	218.999.119	218.999.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.27mm	218.999.127	218.999.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	218.129.090	218.129.090S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.00mm	218.129.100	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	218.129.119	218.129.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.23mm	218.129.123	218.129.123S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	218.129.127	218.129.127S


Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.37mm	218.129.137	218.129.137S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	25mm	0.90mm	218.257.090	218.257.090S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	1.27mm	218.197.127	218.197.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	1.47mm	218.197.147	218.197.147S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	6/12.5mm	1.27mm	218.987.127	218.987.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	6/12.5mm	1.47mm	218.987.147	218.987.147S



Length	Width	Thickness	REF non-sterile	REF sterile
60mm	12.5mm	0.90mm	218.126.090	218.126.090S

Overview references

Length	Width	Thickness				
		0.90mm	1.00mm	1.19mm	1.27mm	1.37mm
<b>Performance Blades</b>						
110mm	25mm				268.251.127 268.251.127S	
110mm	19mm				268.191.127 268.191.127S	
90mm	25mm	268.259.090 268.259.090S	268.259.100 268.259.100S		268.259.127 268.259.127S	
90mm	19mm			268.199.119 268.199.119S	268.199.127 268.199.127S	
90mm	12.5/19mm				268.999.127 268.999.127S	
90mm	12.5mm	268.129.090 268.129.090S	268.129.100 268.129.100S	268.129.119 n/a	268.129.127 268.129.127S	
80mm	12.5mm				268.128.127 268.128.127S	
80mm	8mm				268.088.127 268.088.127S	
70mm	25mm	268.257.090 268.257.090S				268.257.137 n/a
70mm	19mm	268.197.090 268.197.090S				
70mm	12.5mm	268.127.090 268.127.090S			268.127.127 268.127.127S	

Note: sterile version with S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
110mm	25mm	1.27mm	268.251.127	268.251.127S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
110mm	19mm	1.27mm	268.191.127	268.191.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	0.90mm	268.259.090	268.259.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	268.259.100	268.259.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	268.259.127	268.259.127S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	268.199.119	268.199.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	268.199.127	268.199.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.27mm	268.999.127	268.999.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	268.129.090	268.129.090S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.00mm	268.129.100	268.129.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	268.129.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	268.129.127	268.129.127S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	12.5mm	1.27mm	268.128.127	268.128.127S

Performance




Length	Width	Thickness	REF non-sterile	REF sterile
80mm	8mm	1.27mm	268.088.127	268.088.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	25mm	0.90mm	268.257.090	268.257.090S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	25mm	1.37mm	268.257.137	n/a




Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	0.90mm	268.197.090	268.197.090S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12.5mm	0.90mm	268.127.090	268.127.090S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12.5mm	1.27mm	268.127.127	268.127.127S

### Overview references

Length	Width	Thickness						Description	
		0.85mm	1.00mm	1.10mm	1.20mm	1.23mm	1.40mm		
<b>Reciprocating Blades</b>									
85mm	11.6mm				232.785.120 232.785.120S			Reciprocating aggressive	
80mm	10mm						232.980.990 232.980.990S	Reciprocating (file)	
70mm	17mm			232.177.110 232.177.110S		232.177.123 232.177.123S		Reciprocating	
70mm	12mm			232.270.110 232.270.110S				Reciprocating (double-sided)	
60mm	9mm		232.160.100 232.160.100S					Reciprocating	
55mm	10mm	232.755.085 232.755.085S						Reciprocating aggressive	
32mm	6mm			232.132.110 232.132.110S				Reciprocating (sternum)	
cementless		n/a 232.018.275S							Reciprocating Keel Blade
cemented		n/a 232.018.290S							Reciprocating Keel Blade
cementless		n/a 232.000.001S							Reciprocating Keel Blade Set
cemented		n/a 232.000.002S							Reciprocating Keel Blade Set

Note: sterile version with S

### Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	232.785.120	232.785.120S

Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	10mm	1.40mm	232.980.990	232.980.990S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	17mm	1.10mm	232.177.110	232.177.110S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	17mm	1.23mm	232.177.123	232.177.123S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12mm	1.10mm	232.270.110	232.270.110S

Reciprocating




Length	Width	Thickness	REF non-sterile	REF sterile
60mm	9mm	1.00mm	232.160.100	232.160.100S



Length	Width	Thickness	REF non-sterile	REF sterile
55mm	10mm	0.85mm	232.755.085	232.755.085S



Length	Width	Thickness	REF non-sterile	REF sterile
32mm	6mm	1.10mm	232.132.110	232.132.110S



Length	Width	Thickness	REF non-sterile	REF sterile
200mm	6mm	1.10mm	n/a	232.018.275S

cementless

n/a

232.018.275S

**Reciprocating**



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	232.018.290S


**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	218.129.090	218.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	232.785.120	232.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	232.018.275S

\* Keel Blade Set only available sterile

**Keel Blade Set cementless**

**REF non-sterile**

**REF sterile**

n/a

232.000.001S


**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	218.129.090	218.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	232.785.120	232.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	232.018.290S

\* Keel Blade Set only available sterile

**Keel Blade Set cemented**

**REF non-sterile**

**REF sterile**

n/a

232.000.002S

Overview references

Length	Width	Thickness		
		0.60mm	0.80mm	1.20mm
<b>Trauma Blades</b>				
70mm	27mm		208.2770.080 208.2770.080S	208.2770.120 208.2770.120S
70mm	18mm		208.1870.080 208.1870.080S	
60mm	30mm		208.3060.080 208.3060.080S	
48mm	27mm	208.2748.060 208.2748.060S		
48mm	14mm	208.1448.060 208.1448.060S		
35mm	9mm		208.0935.080 208.0935.080S	
25mm	14mm	208.1425.060 208.1425.060S		
25mm	6mm	208.0625.060 208.0625.060S		
24mm	9mm		208.0924.080 208.0924.080S	

Note: sterile version with S

Trauma



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	0.80mm	208.2770.080	208.2770.080S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	1.20mm	208.2770.120	208.2770.120S

Trauma



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	18mm	0.80mm	208.1870.080	208.1870.080S



Length	Width	Thickness	REF non-sterile	REF sterile
60mm	30mm	0.80mm	208.3060.080	208.3060.080S



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	27mm	0.60mm	208.2748.060	208.2748.060S



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	14mm	0.60mm	208.1448.060	208.1448.060S

Trauma



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	9mm	0.80mm	208.0935.080	208.0935.080S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	14mm	0.60mm	208.1425.060	208.1425.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	208.0625.060	208.0625.060S



Length	Width	Thickness	REF non-sterile	REF sterile
24mm	9mm	0.80mm	208.0924.080	208.0924.080S

Overview references

Length	Width	Thickness			
		0.40mm	0.50mm	0.60mm	0.80mm
<b>Small Bone Blades</b>					
45mm	14mm			281.1445.060 281.1445.060S	
45mm	13mm			281.1345.060 281.1345.060S	
41mm	19mm	281.1941.040 281.1941.040S		281.1941.060 281.1941.060S	
40mm	15mm		281.1540.050 281.1540.050S		
39mm	13mm	281.1339.040 281.1339.040S			
35mm	13mm	281.1335.040 281.1335.040S			
35mm	9mm		281.0935.050 281.0935.050S	281.0935.060 281.0935.060S	
31mm	9mm	281.0931.040 281.0931.040S		281.0931.060 281.0931.060S	
31mm	7mm	281.0731.040 281.0731.040S			
31mm	5mm	281.0531.040 281.0531.040S			
25mm	14mm			n/a 281.1425.060S	
25mm	9mm	281.0925.040 281.0925.040S			
25mm	5mm		281.0525.050 281.0525.050S		
21mm	7mm	281.0721.040 281.0721.040S			
18mm	10mm				281.1001.080 n/a
18mm	9mm	281.0918.040 281.0918.040S			
18mm	5mm	281.0518.040 281.0518.040S			
15mm	7mm	281.0715.040 281.0715.040S			

Note: sterile version with S

Small Bone



Length	Width	Thickness	REF non-sterile	REF sterile
45mm	14mm	0.60mm	281.1445.060	281.1445.060S

Small Bone



Length	Width	Thickness	REF non-sterile	REF sterile
45mm	13mm	0.60mm	281.1345.060	281.1345.060S



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	19mm	0.40mm	281.1941.040	281.1941.040S



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	19mm	0.60mm	281.1941.060	281.1941.060S



Length	Width	Thickness	REF non-sterile	REF sterile
40mm	15mm	0.50mm	281.1540.050	281.1540.050S

Small Bone



Length	Width	Thickness	REF non-sterile	REF sterile
39mm	13mm	0.40mm	281.1339.040	281.1339.040S



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	13mm	0.40mm	281.1335.040	281.1335.040S



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	9mm	0.50mm	281.0935.050	281.0935.050S



Length	Width	Thickness	REF non-sterile	REF sterile
35mm	9mm	0.60mm	281.0935.060	281.0935.060S

Small Bone



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	9mm	0.40mm	281.0931.040	281.0931.040S



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	9mm	0.60mm	281.0931.060	281.0931.060S



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	7mm	0.40mm	281.0731.040	281.0731.040S



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	5mm	0.40mm	281.0531.040	281.0531.040S

**Small Bone**



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	14mm	0.60mm	n/a	281.1425.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	9mm	0.40mm	281.0925.040	281.0925.040S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	5mm	0.50mm	281.0525.050	281.0525.050S



Length	Width	Thickness	REF non-sterile	REF sterile
21mm	7mm	0.40mm	281.0721.040	281.0721.040S

**Small Bone**



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	10mm	0.80mm	281.1001.080	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	9mm	0.40mm	281.0918.040	281.0918.040S



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	5mm	0.40mm	281.0518.040	281.0518.040S



Length	Width	Thickness	REF non-sterile	REF sterile
15mm	7mm	0.40mm	281.0715.040	281.0715.040S

# Fits to Synthes® (AO)

## Overview references

Length	Width	Thickness							
Premium Blades		0.90mm	1.00mm	1.18mm	1.19mm	1.23mm	1.27mm	1.37mm	1.47mm
90mm	25mm	215.259.090 215.259.090S	215.259.100 215.259.100S	215.259.118 215.259.118S	215.259.119 215.259.119S	215.259.123 215.259.123S	215.259.127 215.259.127S	215.259.137 215.259.137S	215.259.147 215.259.147S
90mm	22mm						215.229.127 215.229.127S		
90mm	19mm	215.199.090 215.199.090S	215.199.100 215.199.100S		215.199.119 215.199.119S		215.199.127 215.199.127S	215.199.137 215.199.137S	215.199.147 215.199.147S
90mm	12.5/19mm	215.999.090 215.999.090S					215.999.127 215.999.127S		
90mm	12.5mm	215.129.090 215.129.090S	215.129.100 215.129.100S	215.129.118 215.129.118S	215.129.119 215.129.119S	215.129.123 215.129.123S	215.129.127 215.129.127S	215.129.137 215.129.137S	215.129.147 215.129.147S
70mm	25mm	215.257.090 215.257.090S	215.257.100 215.257.100S						
70mm	6/12.5mm						215.987.127 215.987.127S		
60mm	12.5mm	215.126.090 215.126.090S							

Note: sterile version with S

## Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	0.90mm	215.259.090	215.259.090S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	215.259.100	215.259.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.18mm	215.259.118	215.259.118S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.19mm	215.259.119	215.259.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.23mm	215.259.123	215.259.123S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	215.259.127	215.259.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.37mm	215.259.137	215.259.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.47mm	215.259.147	215.259.147S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	22mm	1.27mm	215.229.127	215.229.127S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	0.90mm	215.199.090	215.199.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.00mm	215.199.100	215.199.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	215.199.119	215.199.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	215.199.127	215.199.127S

Premium




Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.37mm	215.199.137	215.199.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.47mm	215.199.147	215.199.147S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	0.90mm	215.999.090	215.999.090S




Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.27mm	215.999.127	215.999.127S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	215.129.090	215.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.00mm	215.129.100	215.129.100S




Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.18mm	215.129.118	215.129.118S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	215.129.119	215.129.119S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.23mm	215.129.123	215.129.123S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	215.129.127	215.129.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.37mm	215.129.137	215.129.137S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.47mm	215.129.147	215.129.147S

Premium



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	25mm	0.90mm	215.257.090	215.257.090S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	25mm	1.00mm	215.257.100	215.257.100S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	6/12.5mm	1.27mm	215.987.127	215.987.127S



Length	Width	Thickness	REF non-sterile	REF sterile
60mm	12.5mm	0.90mm	215.126.090	215.126.090S

Overview references

Length	Width	Thickness				
		0.90mm	1.00mm	1.19mm	1.27mm	1.37mm
<b>Performance Blades</b>						
110mm	19mm				265.191.127 265.191.127S	
90mm	25mm	265.259.090 265.259.090S	265.259.100 265.259.100S	265.259.119 265.259.119S	265.259.127 265.259.127S	265.259.137 265.259.137S
90mm	19mm		265.199.100 265.199.100S	265.199.119 265.199.119S	265.199.127 n/a	
90mm	12.5/19mm				265.999.127 265.999.127S	
90mm	12.5mm	265.129.090 265.129.090S		265.129.119 n/a	265.129.127 265.129.127S	
80mm	12.5mm				265.128.127 265.128.127S	
80mm	8mm				265.088.127 265.088.127S	
70mm	19mm	265.197.090 265.197.090S				

Note: sterile version with S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
110mm	19mm	1.27mm	265.191.127	265.191.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	0.90mm	265.259.090	265.259.090S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.00mm	265.259.100	265.259.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.19mm	265.259.119	265.259.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.27mm	265.259.127	265.259.127S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	25mm	1.37mm	265.259.137	265.259.137S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.00mm	265.199.100	265.199.100S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.19mm	265.199.119	265.199.119S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	19mm	1.27mm	265.199.127	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5/19mm	1.27mm	265.999.127	265.999.127S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	265.129.090	265.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.19mm	265.129.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	1.27mm	265.129.127	265.129.127S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	12.5mm	1.27mm	265.128.127	265.128.127S

Performance



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	8mm	1.27mm	265.088.127	265.088.127S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	19mm	0.90mm	265.197.090	265.197.090S

Overview references

Length	Width	Thickness					Description
		1.00mm	1.10mm	1.20mm	1.30mm	1.40mm	
<b>Reciprocating Blades</b>							
88mm	14mm			245.188.120 245.188.120S			Reciprocating aggressive
85mm	11.6mm			235.785.120 235.785.120S			Reciprocating aggressive
82mm	10mm	245.182.100 245.182.100S					Reciprocating aggressive
80mm	10mm					235.980.990 235.980.990S	Reciprocating (file)
77mm	12mm		235.177.110 235.177.110S				Reciprocating
70mm	12mm		235.270.110 235.270.110S				Reciprocating (double-sided)
70mm	10mm				235.170.130 235.170.130S		Reciprocating
55mm	10mm	245.155.100 n/a					Reciprocating aggressive
cementless		n/a 235.018.275S					Reciprocating Keel Blade
cemented		n/a 235.018.290S					Reciprocating Keel Blade
cemented		n/a 235.016.250S					Reciprocating Keel Blade
cementless		n/a 235.000.001S					Reciprocating Keel Blade Set
cemented		n/a 235.000.002S					Reciprocating Keel Blade Set

Note: sterile version with S

Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
88mm	14mm	1.20mm	245.188.120	245.188.120S

Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	235.785.120	235.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
82mm	10mm	1.00mm	245.182.100	245.182.100S



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	10mm	1.40mm	235.980.990	235.980.990S



Length	Width	Thickness	REF non-sterile	REF sterile
77mm	12mm	1.10mm	235.177.110	235.177.110S

Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	12mm	1.10mm	235.270.110	235.270.110S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	10mm	1.30mm	235.170.130	235.170.130S



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	10mm	1.00mm	245.155.100	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	235.018.275S

Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	235.018.290S



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	235.016.250S

**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	215.129.090	215.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	235.785.120	235.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cementless		n/a	235.018.275S

\* Keel Blade Set only available sterile

**Keel Blade Set cementless**

**REF non-sterile**

**REF sterile**

n/a

235.000.001S

**Keel Blade Set**



Length	Width	Thickness	REF non-sterile	REF sterile
90mm	12.5mm	0.90mm	215.129.090	215.129.090S



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	235.785.120	235.785.120S



Length	Width	Thickness	REF non-sterile	REF sterile
	cemented		n/a	235.018.290S

\* Keel Blade Set only available sterile

**Keel Blade Set cemented**

**REF non-sterile**

**REF sterile**

n/a

235.000.002S

Overview references

Length	Width	Thickness			
		0.60mm	0.80mm	1.00mm	1.20mm
<b>Trauma Blades</b>					
70mm	27mm		206.2770.080 206.2770.080S	206.2770.100 206.2770.100S	206.2770.120 206.2770.120S
70mm	27mm				206.2770.120G n/a
70mm	18mm		206.1870.080 206.1870.080S	206.1870.100 206.1870.100S	206.1870.120 n/a
50mm	20mm		206.2050.080 n/a		
48mm	27mm	206.2748.060 206.2748.060S			
48mm	14mm	206.1448.060 206.1448.060S			
48mm	10mm	206.1048.060 206.1048.060S			
25mm	14mm	206.1425.060 206.1425.060S			
25mm	10mm	206.1025.060 206.1025.060S			
25mm	6mm	206.0625.060 206.0625.060S			

Note: sterile version with S

Trauma



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	0.80mm	206.2770.080	206.2770.080S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	1.00mm	206.2770.100	206.2770.100S


Trauma




Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	1.20mm	206.2770.120	206.2770.120S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	27mm	1.20mm	206.2770.120G	n/a




Length	Width	Thickness	REF non-sterile	REF sterile
70mm	18mm	0.80mm	206.1870.080	206.1870.080S



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	18mm	1.00mm	206.1870.100	206.1870.100S

Trauma



Length	Width	Thickness	REF non-sterile	REF sterile
70mm	18mm	1.20mm	206.1870.120	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	20mm	0.80mm	206.2050.080	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	27mm	0.60mm	206.2748.060	206.2748.060S



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	14mm	0.60mm	206.1448.060	206.1448.060S

Trauma



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	10mm	0.60mm	206.1048.060	206.1048.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	14mm	0.60mm	206.1425.060	206.1425.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	10mm	0.60mm	206.1025.060	206.1025.060S




Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	206.0625.060	206.0625.060S

Overview references

Length	Width	Thickness			
		0.40mm	0.50mm	0.60mm	0.80mm
<b>Small Bone Blades</b>					
50mm	27mm			285.2750.060 285.2750.060S	
50mm	19mm			285.1950.060 285.1950.060S	
50mm	14mm			285.1450.060 285.1450.060S	
50mm	10mm			285.1050.060 285.1050.060S	
48mm	27mm			296.0601.060 n/a	
40mm	15mm		285.1540.050F 285.1540.050FS		
40mm	11mm		285.1140.050F 285.1140.050FS		
31mm	11mm	285.1131.040F 285.1131.040FS			
31mm	7mm	285.0731.040F 285.0731.040FS			
27mm	11mm	285.1127.040F n/a			
25mm	14mm			285.1425.060 285.1425.060S	
25mm	10mm			285.1025.060 285.1025.060S	
25mm	6mm			285.0625.060 285.0625.060S	
22mm	9mm	285.0922.040F 285.0922.040FS			
21mm	13mm	285.1321.040F 285.1321.040FS			
21mm	9mm	285.0921.040F 285.0921.040FS			
18mm	5mm	285.0518.040F 285.0518.040FS			
15mm	15mm	285.1515.040F 285.0518.040FS			
15mm	11mm	285.1115.040F 285.1115.040FS			
15mm	7mm	285.0715.040F 285.0715.040FS			
16mm	10mm				285.1001.080 n/a

Note: sterile version with S

Small Bones Fits to: Colibri, UNIUM™



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	27mm	0.60mm	285.2750.060	285.2750.060S



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	19mm	0.60mm	285.1950.060	285.1950.060S



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	14mm	0.60mm	285.1450.060	285.1450.060S



Length	Width	Thickness	REF non-sterile	REF sterile
50mm	10mm	0.60mm	285.1050.060	285.1050.060S

Small Bones Fits to: Colibri, UNIUM™



Length	Width	Thickness	REF non-sterile	REF sterile
48mm	27mm	0.60mm	295.0601.060	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
40mm	15mm	0.50mm	285.1540.050F	285.1540.050FS



Length	Width	Thickness	REF non-sterile	REF sterile
40mm	11mm	0.50mm	285.1140.050F	285.1140.050FS



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	11mm	0.40mm	285.1131.040F	285.1131.040FS

Small Bones Fits to: Colibri, UNIUM™



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	7mm	0.40mm	285.0731.040F	285.0731.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
27mm	11mm	0.40mm	285.1127.040F	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	14mm	0.60mm	285.1425.060	285.1425.060S



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	10mm	0.60mm	285.1025.060	285.1025.060S

Small Bones Fits to: Colibri, UNIUM™



Length	Width	Thickness	REF non-sterile	REF sterile
25mm	6mm	0.60mm	285.0625.060	285.0625.060S



Length	Width	Thickness	REF non-sterile	REF sterile
22mm	9mm	0.40mm	285.0922.040F	285.0922.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
21mm	13mm	0.40mm	285.1321.040F	285.1321.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
21mm	9mm	0.40mm	285.0921.040F	285.0921.040FS

Small Bones Fits to: Colibri, UNIUM™



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	5mm	0.40mm	285.0518.040F	285.0518.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
15mm	15mm	0.40mm	285.1515.040F	285.1515.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
15mm	11mm	0.40mm	285.1115.040F	285.1115.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
15mm	7mm	0.40mm	285.0715.040F	285.0715.040FS

Small Bones Fits to: Colibri, UNIUM™



Length	Width	Thickness	REF non-sterile	REF sterile
16mm	10mm	0.80mm	285.1001.080	n/a

Overview references

Length	Width	Thickness	
Small Bone Blades		0.40mm	0.60mm
41mm	19mm		286.1941.060 286.1941.060S
31mm	11mm	286.1131.040F 286.1131.040FS	
31mm	7mm	286.0731.040F 286.0731.040FS	
27mm	14mm		286.1427.060 286.1427.060S
27mm	10mm		286.1027.060 286.1027.060S
27mm	6mm		286.0627.060 286.0627.060S
22mm	13mm	286.1322.040F 286.1322.040FS	
22mm	9mm	286.0922.040F 286.0922.040FS	
18mm	5mm	286.0518.040F 286.0518.040FS	
18mm	3mm	286.0318.040F 286.0318.040FS	
15mm	11mm	286.1115.040F 286.1115.040FS	
15mm	7mm	286.0715.040F 286.0715.040FS	

Note: sterile version with S

Small Bones Fits to: Electric Pen Drive



Length	Width	Thickness	REF non-sterile	REF sterile
41mm	19mm	0.60mm	286.1941.060	286.1941.060S

Small Bones Fits to: Electric Pen Drive



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	11mm	0.40mm	286.1131.040F	286.1131.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
31mm	7mm	0.40mm	286.0731.040F	286.0731.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
27mm	14mm	0.60mm	286.1427.060	286.1427.060S



Length	Width	Thickness	REF non-sterile	REF sterile
27mm	10mm	0.60mm	286.1027.060	286.1027.060S

Small Bones Fits to: Electric Pen Drive



Length	Width	Thickness	REF non-sterile	REF sterile
27mm	6mm	0.60mm	286.0627.060	286.0627.060S



Length	Width	Thickness	REF non-sterile	REF sterile
22mm	13cm	0.40mm	286.1322.040F	286.1322.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
22mm	9mm	0.40mm	286.0922.040F	286.0922.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	5mm	0.40mm	286.0518.040F	286.0518.040FS

Small Bones Fits to: Electric Pen Drive



Length	Width	Thickness	REF non-sterile	REF sterile
18mm	3mm	0.40mm	286.0318.040F	286.0318.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
15mm	11mm	0.40mm	286.1115.040F	286.1115.040FS



Length	Width	Thickness	REF non-sterile	REF sterile
15mm	7mm	0.40mm	286.0715.040F	286.0715.040FS

# Fits to Tuke Saw™

Length	Width	Thickness		Description
<b>Reciprocating Blades</b>				
		<b>0.80mm</b>	<b>1.00mm</b>	
87mm	12mm	236.287.080 n/a		Reciprocating (double-sided)
80mm	10mm		236.180.100 n/a	Reciprocating (one-sided)
66mm	16.6mm	236.266.080 236.266.080S		Reciprocating (double-sided)
65mm	12mm	236.265.080 n/a		Reciprocating (double-sided)

Note: sterile version with S

## Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
87mm	12mm	0.80mm	236.287.080	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
80mm	10mm	1.00mm	236.180.100	n/a

Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
66mm	16.6mm	0.80mm	236.266.080	236.266.080S



Length	Width	Thickness	REF non-sterile	REF sterile
65mm	12mm	0.80mm	236.265.080	n/a

# Fits to Zimmer®

## Overview references

Length	Width	Thickness				
Performance Blades		0.90mm	1.00mm	1.19mm	1.27mm	1.37mm
80/100mm	25mm	294.0604.090 n/a	294.0604.100 n/a	294.0604.119 n/a	294.0604.127 294.0604.127S	294.0604.137 n/a
80/100mm	19mm	294.0603.090 n/a	294.0603.100 n/a	294.0603.119 294.0603.119S	294.0603.127 294.0603.127S	294.0603.137 n/a
80/100mm	12.5mm	294.0601.090 n/a	294.0601.100 n/a	294.0601.119 n/a	294.0601.127 294.0601.127S	294.0601.137 294.0601.137S
80/100mm	8mm				294.0602.127 294.0602.127S	

Note: sterile version with S

## Performance Fits to: X-Series



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	25mm	0.90mm	294.0604.090	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	25mm	1.00mm	294.0604.100	n/a

Performance Fits to: X-Series




Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	25mm	1.19mm	294.0604.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	25mm	1.27mm	294.0604.127	294.0604.127S



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	25mm	1.37mm	294.0604.137	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	19mm	0.90mm	294.0603.090	n/a

Performance Fits to: X-Series



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	19mm	1.00mm	294.0603.100	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	19mm	1.19mm	294.0603.119	294.0603.119S



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	19mm	1.27mm	294.0603.127	294.0603.127S



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	19mm	1.37mm	294.0603.137	n/a

Performance Fits to: X-Series



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	12.5mm	0.90mm	294.0601.090	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	12.5mm	1.00mm	294.0601.100	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	12.5mm	1.19mm	294.0601.119	n/a



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	12.5mm	1.27mm	294.0601.127	294.0601.127S

Performance Fits to: X-Series



Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	12.5mm	1.37mm	294.0601.137	294.0601.137S




Length	Width	Thickness	REF non-sterile	REF sterile
80/100mm	8mm	1.27mm	294.0602.127	294.0602.127S

Overview references

Length	Width	Thickness	Description
<b>Reciprocating Blades</b>		<b>1.20mm</b>	
85mm	11.6mm	234.785.120 234.785.120S	Reciprocating aggressive

Note: sterile version with S

Reciprocating



Length	Width	Thickness	REF non-sterile	REF sterile
85mm	11.6mm	1.20mm	234.785.120	234.785.120S

# Fits to Bojin®

## Overview references

Length	Width	Thickness
<b>Premium Blades</b>		<b>1.27mm</b>
100mm	19mm	298.001.127 n/a

Note: sterile version with S

## Premium



Length	Width	Thickness	REF non-sterile	REF sterile
100mm	19mm	1.27mm	298.001.127	n/a

# About Gomina



Gomina is a leading manufacturer of medical products for bone surgery. We develop and produce our surgical instruments exclusively in Switzerland and market them worldwide. Hightech products from the heart of the Swiss mountains.

Gomina has pioneered the development of effective tooth geometry for oscillating saw blades using a combination of sharpened cutting and raker teeth. This unique combination of saw teeth results in outstanding cutting performance. Gomina's tooth geometry and arrangement minimises any lateral swerve during cutting, simplifying handling for the clinician.

Gomina's Keel Blade with its patented, innovative tooth geometry enables the accurate preparation of grooves without bone fragments.

## 100% SWISS MADE

The high requirements of saw technology commanded by surgeons form the basis for our innovation. Today our innovative medical products are used routinely in many knee, hip, foot and hand surgeries to the benefit of patients.

## PRECISION AND SWISS-QUALITY

Since the founding of Gomina in 1974, we have made it our business to be precise.

Since then, Gomina has become a leader in the production of micro-mechanical precision components.

Precision that is used around the world in medical technology and sensor technology.


## PRODUCTION

From the initial blank to the packaged finished unsterile product, Gomina carries out all processing steps in-house, except the sterilization, in compliance with the Swiss Quality Label.

## R&D

Our engineering department converts visionary ideas into innovative products. Thanks to the latest technologies we are able to fulfill special customer requests. In close cooperation with the surgeons we are constantly expanding our product portfolio.

# Certificates



In compliance with the Swiss Quality Label and the standards of the Medical Device Regulation established by the European Union, Gomina fulfills and ensures the strict requirements for safe and efficient medical products. All Gomina products are MDR-certified. Furthermore we fulfill the quality standards within the ISO 13485.

**Our saw blades are certified for reusability**, except for the Reciprocating Keel Blade. The valid certificates, instructions for use, processing and cleaning documents can be found online on [www.gomina.ch](http://www.gomina.ch)

# Compatibility


Gomina is a renowned manufacturer of high quality surgical saw blades for the use in modern orthopaedic, traumatological and reconstructive surgery. Our products are compatible with the power systems of various third party manufacturers established on the market. The compatibility of our replacement saw blades with the different power systems is indicated in this catalogue and on our packaging labels.

Our products are safe and effective to use in, but are not limited to, the power systems of Aesculap, Arthrex, Aygun, Conmed Linvatec Hall, DePuy Synthes, DeSoutter Medical, MatOrtho, MicroAire, Smith & Nephew, Bojin, Stryker and Zimmer Biomet. All trademarks mentioned in this document are registered trademarks of the respective companies. Gomina has no association with or is in any way affiliated with the mentioned companies. Data is correct at the time of publishing and its subject to change without notice.

In the past, uncertainty arose among some users, whether or not it is permitted to use Gomina saw blades (third party product) with their power systems, or if they are obliged to use the „original accessories“ provided by the manufacturer of the power system. This is due to ambiguous statements in the IFUs of some power systems or to misleading information given by the employees of individual power system manufacturers or distributors.

In this context, we would like to confirm that the use of Gomina's surgical saw blades in the power systems specified by us is definitely legally permitted and safe.

For any questions or additional information, please contact us.



# Contact



Gomina AG  
Raiftstrasse 4  
3989 Niederwald, Switzerland  
Phone +41 27 970 12 70  
Fax +41 27 970 12 75  
[www.gomina.ch](http://www.gomina.ch)  
[sales@gomina.ch](mailto:sales@gomina.ch)

**CE** 0482



GOMINA Deutschland GmbH  
Ludwig-Vierling-Strasse 7  
89233 Neu-Ulm, Germany  
[info@gomina.de](mailto:info@gomina.de)

Version 2026-03  
Last modification: 24.03.2026