


# Extraction set for hip and knee revisions

## Universal system



# Handle and extraction bolt


The extraction bolt can be attached on the top or bottom of the handle.



## Usable length of chisels

The handle with a quick clamping mechanism offers three fixation positions for the chisels.

In clamping position 2, the working length of the chisel is extended by 15 mm compared to position 1, and by a further 15 mm in position 3. In total, the working length can be extended by 30 mm from position 1 to position 3.



Usable length =  
L [mm]

# Universal chisels

The design and dimensions of the chisels fit the most common hip stems.  
As single-use instruments, the chisels are supplied in sterile packaging.

| Reference (sterile)   | Article Description          | L [mm] | W [mm] | T [mm] |
|---|------------------------------|--------|--------|--------|
|  <b>401.101S</b>   | Chisel straight              | 143.0  | 6.0    | 1.19   |
|  <b>401.102S</b>   | Chisel straight              | 83.0   | 6.0    | 1.19   |
|  <b>401.201S</b> | Chisel lateral curved        | 116.0  | 5.8    | 1.19   |
|  <b>401.203S</b> | Chisel medial curved         | 114.5  | 5.4    | 1.19   |
|  <b>401.403S</b> | Chisel ventral for right hip | 118.4  | 9.5    | 1.19   |
|  <b>401.404S</b> | Chisel dorsal for right hip  | 118.4  | 9.5    | 1.19   |

# Universal chisels

The design and dimensions of the chisels fit the most common unicondylar knee systems.

As single-use instruments, the chisels are supplied in sterile packaging.

| Reference (sterile)   | Article Description | L [mm] | W [mm] | T [mm] |
|---|---------------------|--------|--------|--------|
|  <b>401.401S</b>   | Chisel straight     | 100.0  | 5.0    | 1.19   |
|  <b>401.402S</b>   | Chisel straight     | 100.0  | 5.0    | 1.19   |
|  <b>401.103S</b> | Chisel straight     | 100.0  | 5.0    | 1.19   |
|  <b>401.104S</b> | Chisel straight     | 100.0  | 10.0   | 1.19   |
|  <b>401.105S</b> | Chisel straight     | 100.0  | 15.0   | 1.19   |
|  <b>401.106S</b> | Chisel straight     | 100.0  | 20.0   | 1.19   |
|  <b>401.107S</b> | Chisel straight     | 100.0  | 20.0   | 1.19   |

# Overview of all references



| Reference                           | Article Description          | L [mm] | W [mm] | T [mm] |
|-------------------------------------|------------------------------|--------|--------|--------|
| 400.001                             | Handle incl. extraction bolt |        |        |        |
| 400.100                             | Extraction bolt (spare part) |        |        |        |
| 401.101S                            | Chisel straight              | 143.0  | 6.0    | 1.19   |
| 401.102S                            | Chisel straight              | 83.0   | 6.0    | 1.19   |
| 401.201S                            | Chisel lateral curved        | 116.0  | 5.8    | 1.19   |
| 401.203S                            | Chisel medial curved         | 114.5  | 5.4    | 1.19   |
| 401.403S                            | Chisel ventral for right hip | 118.4  | 9.5    | 1.19   |
| 401.404S                            | Chisel dorsal for right hip  | 118.4  | 9.5    | 1.19   |
| <i>*suitable for knee revisions</i> |                              |        |        |        |
| 401.401S                            | Chisel straight              | 100.0  | 5.0    | 1.19   |
| 401.402S                            | Chisel straight              | 100.0  | 5.0    | 1.19   |
| 401.103S                            | Chisel straight              | 100.0  | 5.0    | 1.19   |
| 401.104S                            | Chisel straight              | 100.0  | 10.0   | 1.19   |
| 401.105S                            | Chisel straight              | 100.0  | 15.0   | 1.19   |
| 401.106S                            | Chisel straight              | 100.0  | 20.0   | 1.19   |
| 401.107S                            | Chisel straight              | 100.0  | 20.0   | 1.19   |

# Contact



Gomina AG  
Raiftstrasse 4  
3989 Niederwald, Switzerland  
Phone +41 27 970 12 70  
Fax +41 27 970 12 75  
[www.gomina.ch](http://www.gomina.ch)  
[sales@gomina.ch](mailto:sales@gomina.ch)

**CE** 0482

**EU REP**

GOMINA Deutschland GmbH  
Ludwig-Vierling-Strasse 7  
89233 Neu-Ulm, Germany  
[info@gomina.de](mailto:info@gomina.de)

Version 2025-09  
Last modification: 23.09.2025

**GOMINA**  
SWISS QUALITY